

# ANA Glass

## Glass Ionomer

**ANA Glass** is the optimal choice for a modern, strong and high quality Glass Ionomer for multi purpose use.

Thanks to the extremely high compressive strength (300 MPa), perfect condensable structure and handling properties we are tempted to call it "the white amalgam". It can be modelled immediately without sticking to the instruments and have a perfect marginal adaption during the placement.

The fluoride release is high, almost double compared to many of the competitive products on the market, another good reason to use ANA Glass, avoiding secondary caries.

### INDICATIONS

- Restorations of Class I, II, III and V
- Core Build-up
- Base/Liner
- Temporary restorations

### PROPERTIES

- Fluoride release 580 ppm
- Compressive strength 300MPa/30 days



### ASSORTMENT

- Package of 25 capsules in shade **A2** or **A3**
- Cap Applicator

<b>FACT</b>	<b>BENEFIT</b>
Non sticky	Will not pull away from prep or stick to instruments
Packable/Condensable	Easy to handle, place and shape
Adjustable consistency	Fits your exact need of firmness
High compressive strength	Assure durability and wear
High Fluoride release	Highest protection from secondary caries
Designed nozzle/tip	Reach difficult areas more easily



Box 1082, SE-262 21 Ängelholm, Sweden  
Location; Framtidsgatan 1b, SE-262 73 Ängelholm  
Int tel: (+46) 431 44 33 60 - Fax: (+46) 431 44 33 99  
mail@nordiskadental.se - www.nordiskadental.se