

# M2000

## Capsule Mixer

M2000 is our basic model, easy to use but still with highest features and quality. Fully computerized and with a microprocessor verifying the continuously cycles per minute.



### TECHNICAL DATA

<b>Model</b>	M2000
<b>Dimensions</b>	300x220x170
<b>Materials</b>	Metal
<b>Noise</b>	65 dBA
<b>Supply voltage</b>	AC 230 or 110 V
<b>Frequency</b>	50 or 60 Hz
<b>Class MD</b>	I rule I2 in Annex IX 93/42/EEC
<b>Console</b>	IP21
<b>7488 Classification</b>	Type I
<b>Operating Mode</b>	30" ON-180" OFF
<b>Max Working Time</b>	30'
<b>Fuses</b>	1xT2.0A 250V

\* Also available with UK-plug



<b>Easy to use -</b>	Place the cap - set the time - ready to go.
<b>Protective transparent cover -</b>	Reduce noise during high speed mixing. Covers arm and capsule and when opening the lid mixing immediately stops.
<b>Multi-use mixer -</b>	Suitable for all types of dental materials in capsules.
<b>Disinfection friendly -</b>	Safe and simple infection control cleaning. Only one surface cover simplify cleaning with any disinfectant solution.
<b>Reinforced Carbon-Fiber Arms with Universal design -</b>	More durable than metal arms and allows for mixing of all types of encapsulated dental materials.
<b>Vibration absorbers -</b>	The mixer will operate quiet and with minimal vibration, advantages that ensures the stability and keeps the mixer in place.
<b>Speed efficient -</b>	From 0 to 4,200 RPM immediately.
<b>Complete mixing -</b>	Universal mixer for amalgam, glass ionomers and other pre-dosed dental materials in capsules and syringes.



**Nordiska Dental**  
 Box 1082, SE-262 21 Ängelholm, Sweden  
 Int.tel: (+46) 431 44 33 60 Fax:(+46) 431 44 33 99  
 info@nordiskadental.se www.nordiskadental.se

# M4000

## Programmable Capsule Mixer

M4000 is our most advanced mixer, fully computerized and programmable. It gives you the ultimate mix-to-mix consistency, thanks to the micro-processor verifying the continuously cycles per minute.

### TECHNICAL DATA

<b>Model</b>	M4000
<b>Dimensions</b>	300x220x170
<b>Materials</b>	Metal
<b>Noise</b>	65 dBA
<b>Supply voltage</b>	AC 230 or 110V
<b>Frequency</b>	50 or 60 Hz
<b>Class MD</b>	I rule I2 in Annex IX93/42/EEC
<b>Console</b>	IP21
<b>7488 Classification</b>	Type I
<b>Operating Mode</b>	30" ON-180" OFF
<b>Max Working Time</b>	30'
<b>Fuses</b>	1xT2.0A 250V

\* Also available with UK-plug



<b>Easy to use -</b>	A touch control system makes the mixer quick and easy to operate. Loading and removal of the capsule is simplified thanks to the capsule holders design.
<b>Protective transparent cover -</b>	Reduce noise during high speed mixing. Covers arm and capsule and when opening the lid mixing immediately stops.
<b>LED Display -</b>	Exact and clear settings.
<b>Multi-use mixer -</b>	Suitable for all types of dental materials in capsules.
<b>Disinfection friendly -</b>	Safe and simple infection control cleaning. Only one surface cover simplify cleaning with any disinfectant solution.
<b>Reinforced Carbon-Fiber Arms with Universal design -</b>	More durable than metal arms and allows for mixing of all types of encapsulated dental materials.
<b>Vibration absorbers -</b>	The mixer will operate quiet and with minimal vibration, advantages that ensures the stability and keeps the mixer in place.
<b>Speed efficient -</b>	From 0 to 4,200 RPM immediately.
<b>Complete mixing -</b>	Universal mixer for amalgam, glass ionomers and other pre-dosed dental materials in capsules and syringes.
<b>Touch Pad and Memory -</b>	Easy setting and programmable with 10 pre-stored mixing times.



Box 1082, SE-262 21 Ängelholm, Sweden  
Int.tel: (+46) 431 44 33 60 Fax:(+46) 431 44 33 99  
info@nordiskadental.se www.nordiskadental.se