



ANA E-Gel™

37% phosphoric acid

ANA E-Gel™ is well balanced thixotropic gel for precise and safe etching of enamel and dentine.

The consistency and viscosity give precise application and is clearly detectable thanks to the contrast in the color. ANA E-Gel™ is easy to rinse off after treatment and no phosphoric acid is left in the cavity.

INDICATIONS

- Etchant for enamel and dentine
- For direct or indirect restorations

PROPERTIES

- Thixotropic
- Detectable
- Easy to rinse off
- 37% phosphoric acid



ASSORTMENT

-ANA E-Gel™, 2 x 2ml +tips 1332000