



# ANA T-fill™

Temporary filling

**ANA T-fill™** is a self-curing, temporary filling and close-up material for many indications.

It is a ready-to-use-paste and have a very high adhesion to dentine but can easily be removed after its purpose. ANA T-fill™ is gentle to the pulp and have great resistance against chewing pressure for a safe closing of the cavity. The material is eugenol free, not to impede composite fillings

## INDICATIONS

- Temporary close-ups

## PROPERTIES

- Self-curing
- Eugenol free
- High adhesion
- Tight marginal seal



## ASSORTMENT

- ANA T-fill™, white, 2 x 30g      1231000