



# ANA Cal™

## Calcium Hydroxide

ANA Cal™ gives a high concentration of hydroxyl ions, thanks to the high content of calcium hydroxide (>38%), still with a smooth consistency. These generate the high pH value of 12.4 which has a pronounced bactericidal effect, i.e. eliminate bacteria. ANA Cal™ is radio-opaque.

High concentration of calcium hydroxide also means long-lasting treatment due to that the paste can release calcium ions for a longer period of time. ANA Cal™ is a ready-to-use product, in air tight syringes, for direct application thru the included needles.

### INDICATIONS

- Temporary root-filling
- Pulp capping
- Pulp protection
- Insulation in deep cavities
- Stepwise excavation



### ASSORTMENT

- Package of 4 syringes, each containing 1,5ml and 12 luer-lock needles

Art.Nr. 123 01 00

FACT	BENEFIT
High concentration, >38%	Powerful, long-lasting effect
Pure calcium hydroxide	Optimal ion release
Smooth consistency	Precise and deep access
Pre filled syringe	Direct, quick and easy application
Multi-purpose	One product - many indications