

A reliable friend

ANA Amalgams






”It isn’t more silver that
makes a superior amalgam.
It’s more copper.”

ANA 2000®

A complete high-quality range of amalgams

The range offers admixed alloy of highest strength, dispersed alloy of smoothness and top clinical handling properties.

Our amalgams are superior in many quality aspects. They have well balanced expansion and high strength. This combination gives tight margins, long lasting fillings and few marginal breakdowns.



Combination of spherical and lathe cut particles, ANA 2000

COMPOSITION

Ag 43,1%
Sn 30,8%
Cu 26,1%

PHYSICAL PROPERTIES

Compressive strength:
1 h 140 MPa
24 h 450 MPa

Static creep: 0.1%

Dimensional change:
24 h +0.02%

ANA 2000[®]
High Copper Admixed Alloy
Non Gamma 2

AVAILABLE IN:

DUETT
CAPSULES
POWDER

ANA 2000

■ FACT	● BENEFIT
High copper	Long-lasting, resistant and strong
Controlled expansion	Excellent marginal adaptation
Lathe cut particles	Crunchy feeling when condensing
Short setting time	Time saving

Duett® – the safety and convenience of capsules without the cost!

The Duett capsule system is environmentally friendly when comparing with other capsule presentations. The reusable capsules are specially designed to be tight, ensuring safety and hygiene. Duett is also cost saving, with reduced cost when comparing with disposable capsules. The convenience of having one box, instead of one for each spill, is indisputable. Dentists conclude that the best homogeneously mix is accomplished with Duett.

DUETT® – UNIQUE PACKAGING

FACT	BENEFIT
Cost effectiveness	Cost saving compared to other types of capsules
Environmentally friendly	Less waste material
Convenience	Easy to dispense and mix
Safety	No loose mercury
Hygiene	Low mercury exposure



PACKAGING SIZES:



Duett® 400
400 tablets and 400 mercury sachets, 8 reusable capsules




Duett® 200
200 tablets and 200 mercury sachets, 4 reusable capsules



Duett® 50
50 tablets and 50 mercury sachets, 2 reusable capsules



NOTE
The alloy tablet contains 390 mg and the mercury sachet 390 mg.



REUSABLE CAPSULE

EASY TO USE

Use the tweezers when charging the capsules with one or two tablets and one or two sachets of mercury. This will not take long and when dosing the capsules in advance you will always have capsules ready to use. The quantity is always the same.

Powder – flexible mixing opportunities

Nordiska alloy is well renowned and therefore speaks for itself. Generations of dentists all over the world have made long lasting fillings of ANA 2000.

POWDER

BENEFIT
Flexible mixing ratio
Price



500 g.



250 g.



30 g.

POWDER

Powder is available in sealed bottles containing 500 g, 250 g or 30 g.

Pre-dosed disposable Capsules

With a precise dosage in each disposable capsule the consistency is guaranteed – consistently!
Today when time is of the essence it is important to have a dispensing system that allows quick and easy handling, our capsules with low mercury exposure makes this possible. Being self activating makes our capsules very convenient.

ANA 2000 -

The ANA 2000 disposable capsules have been optimized reaching yet a higher level of quality. The new special design gives a homogenous mix.

CAPSULES

FACT	BENEFIT
New design	Optimal amalgamation Safe and tight
Pre-dosed	Exact consistency Time saving
Self activating	Easy to mix
Colour coded	Easily differentiated



ANA 2000 CAPSULES

- No. 1 – Green capsule
(400 mg alloy+400 mg mercury)
- No. 2 – Blue capsule
(600 mg alloy+600 mg mercury)
- No. 3 – Grey capsule
(800 mg alloy+800 mg mercury)





Nordiska Dental AB, Box 1082, SE-262 21 Ängelholm, *Location*; Framtidsgatan 1b, SE-262 73 Ängelholm, Sweden,
Tel: +46 431 44 33 60, Fax: +46 431 44 33 99
www.nordiskadental.se, mail@nordiskadental.se